

LARYNGOSCOPES

ARMO diagnostics

Our Armo brand instruments were created in response to medical professionals need for dependable instrumentation in a cost-sensitive environment. Armo laryngoscopes are crafted from the finest quality materials and made to international specifications in state of the art accredited manufacturing facilities.

Laryngoscopes (Standard)



Laryngoscope features

- Stainless steel (416L)
- Laryngoscope Sets are supplied with blades and handle in protective hard case
- McIntosh Set includes McIntosh blades 1 to 4 and a 32mm handle
- Miller Set includes Miller blades 0 to 3 and a 32mm handle
- Sets are supplied excluding batteries

Handle features

- Handles are interchangeable with McIntosh and Miller blades (all sizes)
- 19mm and 32mm diameter handles with sure grip surfacing
- Compatible with most standard hook-on blades that conform to ASTM(F965)

Blade features

- Surgical grade stainless steel
- Autoclavable to 134 °C after removal of bulb
- Compatible with most standard hook-on handles that conform to ASTM(F965)

Ordering Information

Description	GVP Code
McIntosh Set	GVPBD-A0204
Miller Set	GVPBD-A0614
Handle Ø 32mm	GVPBD-A0192
Handle Ø 19mm	GVPBD-A0195
McIntosh Blade #0	GVPBD-A0191
McIntosh Blade #1	GVPBD-A0187
McIntosh Blade #2	GVPBD-A0188
McIntosh Blade #3	GVPBD-A0189
McIntosh Blade #4	GVPBD-A0190
Miller Blade #0	GVPBD-A0615
Miller Blade #1	GVPBD-A0616
Miller Blade #2	GVPBD-A0617
Miller Blade #3	GVPBD-A0618

LARYNGOSCOPES

Fibre Optic (F.O.) Laryngoscope

Laryngoscope features

- Laryngoscope Sets are supplied with blades and handle in protective hard case
- McIntosh Set includes McIntosh blades 1 to 4 and a 32mm handle
- Miller Set includes Miller blades 0 to 3 and a 32mm handle
- Sets and handles are supplied without batteries

Handle features

- Handles are interchangeable with McIntosh and Miller blades (all sizes)
- 32mm diameter handle with sure grip surfacing
- Handles are compatible with most F.O. hook-on blades conforming to the ISO 7376 (green standard)

Blade features

- Surgical grade stainless steel
- Autoclavable to 134 °C
- Handles are compatible with most F.O. hook-on blades conforming to the ISO 7376 (green standard)

Ordering Information

Description	Code
F.O. McIntosh Set	GVPBD-A0201
F.O. Miller Set	GVPBD-A0202
F.O. Handle	GVPBD-A0203
F.O. McIntosh Blade #1	GVPBD-A0716
F.O. McIntosh Blade #2	GVPBD-A0717
F.O. McIntosh Blade #3	GVPBD-A0718
F.O. McIntosh Blade #4	GVPBD-A0719
F.O. Miller Blade #0	GVPBD-A0815
F.O. Miller Blade #1	GVPBD-A0816
F.O. Miller Blade #2	GVPBD-A0817
F.O. Miller Blade #3	GVPBD-A0818
Spare light bulb for F.O.	GVPBD-A0226



Laryngoscope Blade Lengths

Type & Size	Length
McIntosh Blade #0	85mm
McIntosh Blade #1	92mm
McIntosh Blade #2	110mm
McIntosh Blade #3	130mm
McIntosh Blade #4	155mm
Miller Blade #0	78mm
Miller Blade #1	100mm
Miller Blade #2	155mm
Miller Blade #3	195mm